



# Spirobank II<sup>®</sup> Advanced

Advanced handheld diagnostic Spirometer with interpretation.  
Measures 45+ parameters with Oximetry optional. PC connectivity via USB and Bluetooth<sup>®</sup>.

The Spirobank II<sup>®</sup> Advanced is a handheld, portable diagnostic spirometer that offers a comprehensive range of spirometric parameters for more in-depth spirometry testing.

Its pocket size design with rechargeable long-life battery makes it ideal for use in situations where remote testing is required. The user interface is fast and intuitive to operate with an icon-based navigation menu on a backlit, high resolution display. It can be used as a "stand-alone" with onboard memory that is capable of storing up to 10,000 test results. It can also connect wirelessly to a PC (or directly via USB) allowing for real-time spirometry testing, and downloading of test results stored on the device using the free MIR WinspiroPRO<sup>®</sup> PC Software.

The Spirobank II<sup>®</sup> Advanced measures over 45 selectable parameters and gives a PRE and POST bronchodilator comparison. Results are previewed on the display following a test, with a flow/volume graph in addition to the patient's actual spirometry readings for immediate assessment. It also features a traffic light health indicator for easy interpretation of results.

## Re-usable and disposable consumable options

The Spirobank II<sup>®</sup> Advanced can be used with both re-usable and single-use disposable turbines. Re-usable turbines are used with standard adult one-way valve mouthpieces. The single-use disposable turbines have a pre-connected mouthpiece and are individually factory calibrated and packaged.

## KEY FEATURES

- **Easy to use handheld spirometer**  
Navigate easily through the icon-based menus on a clear, high resolution screen
- **Over 45 selectable parameters**  
Advanced handheld diagnostic spirometry ideal for use in any clinical setting
- **Traffic light health indicator**  
A visual way to easily interpret test results
- **Large memory storage**  
Store up to 10,000 tests
- **Bluetooth<sup>®</sup> and USB connectivity to PC**  
Download results stored on the device and conduct real-time spirometry using WinspiroPRO<sup>®</sup> PC software.
- **Optional Pulse Oximetry**



# MIR Spirobank II® Advanced

TECHNICAL SPECIFICATIONS	
Device category	Diagnostic spirometer
Spirometry parameters	FVC, FEV1, FEV1%, FEV1/FVC, FEV1/VC, PEF, FEF25, FEF50, FEF75, FEF25-75, FEF75-85, Lung Age, Extrap. Volume, FET, Time to PEF, FEV0.5, FEV0.5/FVC, FEV0.75, FEV0.75/FVC, FEV2, FEV2/FVC, FEV3, FEV3/FVC, FEV6, FEV1/FEV6, FIVC, FIV1, FIV1/FIVC, PIF, FIF25, FIF50, FIF75, FEF50/FIF50, VC, IVC, IC, ERV, IRV, Rf, VE, VT, tl, tE, VT/tl, tE/tTOT, MVV (measured), MVV (calculated)
Flow measurement	± 16 L/s
Volume accuracy	± 3% or 50 mL
Flow accuracy	± 5% or 200 mL/s
Temperature sensor	Semiconductor (0-45°C), automatic BTPS conversion
Data transmission	USB and Bluetooth® 2.1
Memory storage	Approx. 10,000 tests
Power	3.7 V, 110 mA lithium battery - USB rechargeable
Display	160 x 80 pixel with 6-key membrane keyboard
Dimensions	160x55x25 mm
Weight	145 g (with battery)
<b>Pulse Oximetry Specifications (Upgrade to be purchased separately)</b>	
SpO2 range/accuracy	0-99%SpO2 / ±2% between 70-99%SpO2
Pulse range/accuracy	18-300BPM / ±2BPM or 2% (Whichever is greater)
Oximetry alarms	SpO2 and pulse rate Min-Max, finger/sensor inserted, battery out of energy
Oximetry tests	SpO2. Min/Max/Avg, pulse rate, test duration, T90% [SpO2 time ≤89%], T89% [SpO2 time ≤88%], total SpO2 events, T40 (bradycardia duration with pulse rate <40 BPM). T120 (Tachycardia duration with pulse rate >120 BPM)

## Ordering Information

PART NUMBER	DESCRIPTION
MR-911020E1	MIR Spirobank II® Advanced with re-usable turbine
MR-911020E2	MIR Spirobank II® Advanced with 100 x disposable turbine
MR-910002	Replacement re-usable turbine
MR-910004	Disposable turbines (Box of 50)
881002	One-way valve adult mouthpieces for re-usable turbine (Box of 200)

*For Pulse Oximetry upgrade, please contact the provider*

© 1997-2018. Intermedical (UK) Limited. All rights reserved.

All information, specifications and illustrations contained in this document are based on the latest product information available at the time of publication. Errors and omissions excepted. Intermedical (UK) Limited reserves the right to change specifications and design at any time without prior notice.



Exclusive UK distributor

