



SpiroCal-3000

3 Litre Spirometry Calibration Syringe.
Low cost, accurate calibration syringe for your spirometer

3- Litre calibration syringe suitable for use with:

- ndd EasyOne Spirometer
- ndd EasyOn-PC Spirometer
- MIR Spirometer Range
- Micro Medical Range
- Vitalograph Range

TECHNICAL DATA

Nominal volume	3000 ml
Accuracy	±15 ml
Lubricant	Inteflon, Fin Food Lube
Size	520 x 120 x 110 mm
Weight	2.2kg
Temperature	Operating: 15+40°C Storage: 0+50°C
Humidity	10%+95%

ORDERING INFORMATION

1690000	SpiroCal-3000 3L Syringe
---------	--------------------------

© 1997-2018. Intermedical (UK) Limited. All rights reserved.

All information, specifications and illustrations contained in this document are based on the latest product information available at the time of publication. Errors and omissions excepted. Intermedical (UK) Limited reserves the right to change specifications and design at any time without prior notice.

