

INOGENONE® G3

Portable Oxygen Concentrator

The INOGENONE® G3 is a popular choice amongst oxygen users for its ease of use, stylish and non-medical feel.

With up to 5 pulse flow settings available, the INOGENONE® G3 strikes the perfect balance between those wanting a lightweight solution and the flexibility of delivering oxygen at a more demanding level.



Pulse Flow

1 to 5
Pulse flow settings



Weight

From
4.9lbs
(2.2kgs)



Battery life

Up to
8
Hours
on setting 2*



intermedical
H O M E C A R E

www.portableoxygen.co.uk | 0800 028 2194

Trusted in Respiratory Health

Smart, stylish design with the user in mind.

This compact oxygen solution weighs from only 4.9lbs (2.2kgs) and can be worn easily over the shoulder in its protective black bag. A backpack accessory is also available separately for even greater comfort and independence.

The G3 is powered by either an 8 or 16 cell battery that slots neatly onto the base of the machine with an approximate run-time of 4 or 8 hours respectively (on setting 2).

The user-interface is clear and easy to operate with a white back-lit LCD display that shows the remaining battery operation time in hours/minutes and the current setting.




◀ Backpack accessory is also available for the G3.



◀ The INOGENONE® G3 has two battery options: 8 cell & 16 cell. The battery re-charges when installed and the machine is connected to mains or DC power. A desktop charger is also available to charge any supplementary batteries. (optional accessory)

8 Cell Battery Package

- INOGENONE® G3 main unit
- 8 cell rechargeable battery
- Carry bag with shoulder straps
- AC mains power supply unit
- DC power supply unit
- User manual
- Quick start guide
- Nasal cannula



21.0 (h) x 22.2 (w) x 7.6 (d)
in cm with 8 cell battery installed

16 Cell Battery Package

As above including:

- 16 cell rechargeable battery instead of 8 cell battery

Approximate battery run-time (in hours:minutes)

	Standard Battery	Double Battery
Setting 1	4:30	9:00
Setting 2	4:00	8:00
Setting 3	3:00	6:00
Setting 4	2:00	4:00
Setting 5	1:40	3:20

*Run-time dependent on breathing rate

Concentration	90% -3% /+ 6%
Weight	4.9lbs (2.2kg) - 8 cell batt. 5.7lbs (2.6kg) - 16 cell batt.

Power AC: 100-240v DC: 13.5-15.5v

Noise 39db(A) on setting 2

Batt. charge time 8 cell: 4 hours
16 cell: 8 hours

Warranty 3 years on unit
1 year on batteries, sieve beds & accessories

✈ Airline approved Yes, subject to airline's terms & conditions

~~CPAP~~
~~BiPAP~~ This device may be suitable for use whilst sleeping. This device cannot be used with a CPAP or BiPAP device.