

To place your order please call
0800 028 2194

Or buy online at www.intermedicaldirect.co.uk/homecare

Intermedical (UK) Limited, Homecare Division
Unit 6 Mill Hall Business Estate, Aylesford, Kent, ME20 7JZ, United Kingdom

Trusted in **Respiratory Health**



INOGENONE[®]G4

Optional accessories
for your INOGENONE[®] G4 Portable Oxygen Concentrator

intermedical
H O M E C A R E

INOGENONE® G4 | Accessories List



Highly Recommended

4 cell rechargeable battery

Up to 2:40hrs run-time. Consider purchasing extra battery packs to extend your machines operation time

944010



Highly Recommended

8 cell rechargeable battery

Up to 5 hrs run-time. Consider purchasing extra battery packs to extend your machines operation time

944011



Highly Recommended

Desktop battery charger

Highly recommended if you have more than one battery pack

944013



Carry strap

Replacement shoulder strap that clips on to the INOGENONE® G4 unit or carry bag

944017



Particle filter

Filter housed on the front of the machine recommended to be changed when soiled

944015



AC mains power supply unit

Replacement power supply unit that connects the machine to the mains power supply. Supplied with UK plug

944012



DC 12v power supply unit

Replacement power supply unit that connects the machine to a 12v socket

943015



Carry bag for INOGENONE® G4

Replacement carry bag

944016