

Choosing the right Smokerlyzer



The stop smoking advisor's guide to choosing the best CO monitor.

Need help?

Speak to one of our specialists today on 01732 522444.



pico™

Ideal for use for general smokers in smoking cessation programmes and research.



pico baby™

Ideal for use within the midwifery setting as well as for use with general smokers in smoking cessation programmes and research.



Micro+™

Advanced built-in functionality making it ideal for all round use in breath and ambient air CO monitoring.

Measures ppm and %COHb	✓	✓	✓
For use with pregnant women (measures %FCOhb)	✗	✓	✓
Touchscreen control	✓	✓	✓
No calibration required	✓	✓	✓
5-year warranty	✓	✓	✓
Input of patient data on device	✗	✗	✓
Memory storage	✗	✗	✓
Ambient CO monitoring	✗	✗	✓
Adjustable interpretation ranges	✗	✗	✓
Exports data to PC using free COdata+ Software	✓	✓	✓
PC real-time testing with Smoking in Pregnancy Software	✗	✗	✓
Concentration range	0-150ppm	0-150ppm	0-500ppm

Order Code: 1420004N

**NHS Supply Chain
FDD4407**

Order Code: 1420004NB

**NHS Supply Chain
FDD3223**

Order Code: 1420003N

**NHS Supply Chain
FDD4408**

intermedical
CARDIO RESPIRATORY

Buy on **NHS Supply Chain** or online at www.comonitors.com