


# The simple solution to bring **24hr Holter ECG** services into your GP practice.

All-in-one equipment & remote cardiac  
analysis service from Intermedical.



**intermedical**  
CARDIO RESPIRATORY

A professional Holter ECG system together with interpretation software is notoriously expensive.

The level of expertise and clinical time required to analyse the data, puts offering a holter ECG to patients in a primary care setting out of reach for most practices.



**Our service is designed to make Holter ECG more practical, more affordable and more accessible to primary care - so you can reduce unnecessary referrals to hospital and give patients a faster route to treatment.**

# How it works.

Our service uses the Mortara H3+ Holter ECG device. The PC software that's included facilitates the set-up of the device, offloading and transfer of data to a secure cloud where it is analysed by our partner Express Diagnostics, who are SCST-accredited data analysts. A report is prepared and sent to the practice via secure NHS email.



Patient visits GP with a possible cardiac complaint



Patient is fitted with the Mortara H3+ Holter ECG device in the practice to take home



Patient returns 24 hours later and the data is uploaded for analysis via secure data transfer software



Data is analysed by a cardio physiologist. A report is prepared and returned back to the GP via NHS secure email

The only involvement of the health professional is the correct fitting of the holter and the uploading of data from the device to the cloud where it will be carefully analysed by cardio physiologist.

# How it benefits.

Practical. Cost effective.  
Time-saving. **Life-saving.**

GP surgeries and CCG's employing this practice have the potential to make huge cost and time savings, whilst keeping patients out of hospital with a quality diagnosis and a faster route to treatment.

## For GP Surgeries...

- Affordable access to capital equipment
- No lengthy analysis/interpretation necessary
- Quality assured diagnosis
- Basic clinical expertise needed

## For the CCG...

- Save money by dramatically reducing secondary care referrals
- Reduce wait times & improve patient scoring

## For the Patient...

- Keep patients out of hospital
- Save lives with a faster route to diagnosis



# The technology.

Weighing in at just one ounce (28 grams), the Mortara H3+ Holter ECG is small and lightweight, making it very comfortable for patients. The device is fitted with one AAA battery that is capable of recording non-stop for up to 2 days.

3

channel ECG

28g

in weight

64mm

mm in width

1

year warranty

Unobtrusive miniature design

3-channel ECG recording with pacemaker spike detection

Graphic display to preview ECG waveforms

Only one AAA battery required



 Mortara H3+



Trusted in **Respiratory Health.** Since 1997.

t 01732 522444 e [sales@intermedical.co.uk](mailto:sales@intermedical.co.uk) w [intermedical.co.uk](http://intermedical.co.uk)

In partnership with



Intermedical (UK) Limited | Unit 6 Mill Hall Business Estate | Aylesford | Kent | ME20 7JZ | United Kingdom  
Tel: 01732 522444 Fax: 01732 872883 Email: [sales@intermedical.co.uk](mailto:sales@intermedical.co.uk)

Registered in England & Wales. Company No.3456073

Doc 0420iss01